


### Endoscopic system



Portable endoscopic monoblock with lighting system, IPS screen, insufflator and digital block for working with flexible and rigid endoscopes

Rating: Not Rated Yet

[Ask a question about this product](#)

#### Description

Portable endoscopic monoblock with lighting system, IPS screen, insufflator and digital block for working with flexible and rigid endoscopes

#### Features of the endoscopic monoblock:

- IPS LCD screen for high-quality contrast resolution
- The optional 15-inch screen allows the device to be maximally adapted to clinical needs
- Carry out examinations with rigid and flexible endoscopes, thanks to the adapter
- It is possible to use endoscopes from other manufacturers
- Improved imaging technology and color contrast

#### Specifications

- 12 inch IPS monitor
- Digital matrix CMOS - 1024x720 p
- Auto white balance
- Digital Signal - DVI, HDMI
- Color adjustment - 11 levels
- Saturation adjustment - 11 levels
- Brightness adjustment - 11 levels
- Capturing an image
- Video recording
- Saving images and videos to a memory card
- Data management system
- LED lighting system
- Lighting power 40 W
- Average life of the lighting system 35,000 hours
- Color temperature 5300 K
- Manual dimming
- Air insufflation and water pumping
- Power supply 100 - 240 VAC
- Dimensions (WxHxD) - 375x320x130 mm
- Weight with monitor 9.5 kg

