

Seam and ligature



The needle holder is an endoscopic instrument for grasping the needle when suturing tissues during minimally invasive interventions on the abdominal organs.

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

The needle holder is an endoscopic instrument for grasping the needle when suturing tissues during minimally invasive interventions on the abdominal organs.

NO.	Name	Characteristics
NH0003L	Straight ring handle needle holder	Diameter 5 mm, working length 100 mm
NH0004L	Straight V-handle needle holder	Diameter 5 mm, working length 100 mm
NH0005L	Curved V-Handle Needle Holder	Diameter 5 mm, working length 100 mm
NH0008L	Curved needle holder for ring grip	Diameter 5 mm, working length 100 mm
NH0006L	Needle holder for straight handle	Diameter 5 mm, working length 100 mm
NH0007L	V-handle needle holder	Diameter 5 mm, working length 100 mm