


Cystoscopes



Laparoscopic scissors are an endoscopic instrument with two-edged counter-cutting elements. A knife or laparoscopic scalpel is a classic surgical instrument for tissue dissection.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

The Brightfield cystoscope is a medical device equipped with a sophisticated optical and illumination system with an access channel for microsurgical instruments. The cystoscope is intended for examining the walls of the urethra and bladder, catheterization of the ureters, washing the bladder, conducting minimally invasive surgical procedures (biopsy, coagulation, extraction of foreign bodies.).

Features of Brightfield cystoscopes

- The body of the tube and its constituent elements are made of high grade stainless steel
- The endoscope uses elements of optics from a German manufacturer.
- The optical axes are created using the patented lens shaft technology, which ensures the visibility of the entire field of view.
- The cover and sapphire lens are reliable in operation and never wear out.
- The valve is made of stainless steel, easy to maintain, protected from damage.
- The connectors are compatible with endoscopic equipment from leading manufacturers, including Storz products.
- In urology new methods and technologies, which use the advantages of laser over traditional treatment, continue to appear.

Brightfield Cystoscope Sheath Models

NO.	Field of view	Diameter	Working length
CYS01	0 degrees	4 mm	302 mm
CYS02	30 degrees	4 mm	302 mm
CYS03	70 degrees	4 mm	302 mm
CYS04	0 degrees	2.7 mm	190 mm
CYS05	30 degrees	2.7 mm	190 mm