

## Cystolithotriptor



Cystolithotriptor is an instrument for fragmentation of bladder stones under visual control using mechanical action

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

Cystolithotriptor is an instrument for fragmentation of bladder stones under visual control using mechanical action.

### Features of the Brightfield cystolithotriptor

- The endoscope is made of high quality stainless steel
- Device optical elements of the German optical company.
- The optical axes are created using the patented lens shaft technology, which ensures the visibility of the entire field of view.
- The cover and sapphire lens are reliable in operation and never wear out.
- The valve is made of stainless steel, easy to maintain, protected from damage.
- The connectors are compatible with endoscopic equipment from leading manufacturers, including Storz products.
- Sterilization can take place at low temperatures.

### Brightfield Cystolithotriptor Models

NO.	Field of view	Diameter	Working length
LITH_01.	12 degrees	4 mm	302 mm
LITH_02.	70 degrees, 90 degree viewing angle	4 mm	302 mm
LITH_03.	70 degrees	3 mm	302 mm