

Lighting system U9125 YG-250



U9125 YG-250 Endoscopy Lighting System

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

U9125 YG-250 Endoscopy Lighting System

The light source in the U9125 YG-250 is two halogen lamps. The presence of an additional lamp guarantees a problem-free endoscopic examination or surgery. The modern electrical component of the device provides stable light output indicators.

Main technical characteristics of U9125 YG-250

- Light source - halogen lamp
- Halogen lamp power 250 W
- Voltage 220 V
- Illumination over 1,800,000 lx
- Color temperature 3200 K