


Disconnecting instruments for phonosurgery



Cutting endoscopic instruments for phonosurgery. The universal handle, to which the working tool is attached, makes it possible to set different working lengths.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Cutting endoscopic instruments for phonosurgery. The universal handle, to which the working tool is attached, makes it possible to set different working lengths.

| NO. | Name | Characteristics |
|---------|--|------------------------------|
| KLOOO1L | Lancet knife, straight | Diameter 2 mm, length 250 mm |
| NLOOO1L | Needle, straight, pointed | Diameter 2 mm, length 250 mm |
| HLOOO1L | Hook, tracheoscopic | Diameter 2 mm, length 250 mm |
| HO0000L | Handle for use with tools: KLOOO1, NLOOO1, HLOOO1L | |