

## Otoscope LTI-304



High quality Heine type plastic mini otoscope

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

High quality Heine type plastic mini otoscope

High quality Heine type plastic mini otoscope with extra bright light 2.5v LED illumination bulb illumination provided through fiber optic system. High quality swiveling viewing window with 3x magnification lens, out put for pneumatic membrane test, easily inserts battery facility in handle.

**Supply:** 10 single use mini speculum (2.5 and 4mm) with storage bottle packing in black pouch.