

Otoscope LTI-2010



Pocket instrument for endoscopic diagnostic examinations of the external auditory canal and tympanic membrane with the possibility of pneumotoscopy

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Pocket instrument for endoscopic diagnostic examinations of the external auditory canal and tympanic membrane with the possibility of pneumotoscopy.

The device is compact and easy to use. The handle body and optical head are made of high quality durable plastic. The switch is located on the clip. When the tool is placed in the pocket, it is switched off automatically

The otoscope is equipped with an innovative fiber-optic illumination system that is 2 times superior to standard technologies. The light source is an LED. The power source is 2 AA alkaline batteries (1.5 V). The user can use batteries with similar dimensions and specifications.

The optical system of the otoscope provides standard 3x magnification. The resulting image is displayed on a 4-angle screen. The basic set includes disposable funnels.