

## Laryngoscope LTI-6006



Fiber optic laryngoscope

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

Fiber optic laryngoscope

hi-light fiber optic bundle with 3.5mm diameter for improved light transmission and longer life.

- LED Technology: For up to 40% more light than conventional halogen bulb
- smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing.
- these blades need no maintenance
- both ends of the fiber bundle are protected
- 5V battery power for even more light.
- compatible with all handles to the ISO 7376 (green Standard)
- fits all standard EO. handles (green standard)