

## Laryngoscope LTI-7001



Innovative fiber optic technology for advanced performance

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

Innovative fiber optic technology for advanced performance

- hi-light Fiber Optic bundle with 4.3mm dia. And up to 6500 individual micro-fiber for improved light transmission and longer life.
- LED Technology: For up to 40% more light than conventional halogen bulb
- smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing. :
- these blades need no maintenance
- ever both ends of the fiber bundle are compatible with all handles to the ISO 7376 (green Standard)
- fits all standard EO. handles (green standard)