

## Facial bed BR-273B



4 motor facial bed

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

#### **4 motor facial bed**

Facial bed BR-273B is an extremely comfortable chair/bed. It features 4 motors that control the height, the backrest and the inclination of the seat. The headrest comes with breather hole incorporated and is removable. The armrests are reversable as well as removable. The leg rest is controlled manually.

### Features:

- Sturdy frame
- White PU upholstery
- Headrest breather hole
- Adjustable angle headrest
- Backrest/ footrest can be adjusted to 80 degree
- Height can be adjusted from 55 cm-75 cm by removing controller, it can be tilted to 30 degree
- Bed surface is divided into 4-piece couch