

## Light Source - POWER LED 180



Rating: Not Rated Yet

[Ask a question about this product](#)

Description

REF.-NO.: 103.0180

### Features:

- **NEW:** Optimized colour temperature for natural light (daylight)
- „State of the Art“ due to the latest LED technology
- Light intensity and color temperature comparable with Xenon light source 180 Watts
- No further purchase of overpriced Xenon lamps
- Diode for extreme long life!

### Technical Data:

- **Dim.:** 330 x 310 x 90, Weight: 5 kg
- **Lifetime LED:** > 50.000 hours
- **Power consumption:** 70 VA
- **Power supply:** 115-230 VAC (90-263 V)
- **Frequency:** 50/60 Hz
- **Color temperature:** 5.500 K!
- **Colour Rendering Index:** CRI 95!
- **Safety class:** BF