

Brightfield Healthcare 9000D



15" patient monitor

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

15" patient monitor

This monitor is suitable for adult, pediatric and neonate patient, may monitor the life parameters, such as the electrocardiograph (ECG), the noninvasive blood pressure (NIBP), the oxygen saturation (SpO₂), the respiration rate (Resp), the body temperature (Temp), the invasive blood pressure (IBP), carbon dioxide (CO₂) and so on, can in the same screen demonstrate maximum 11 waveforms and all information of the parameters monitored.

Below shows the monitoring functions of this monitor:

- electrocardiograph (ECG), including: heart rate, 6 channels ECG waveforms, ST segment analysis (can be partially optional and purchased), arrhythmia analysis (can be partially optional and purchased);
- oxygen saturation (SpO₂), including: oxygen saturation, pulse rate, pulse wave;
- noninvasive blood pressure (NIBP), including: Systolic pressure, diastolic pressure, mean pressure;
- body temperature (Temp): 2 channels body temperature data;
- respiration (Resp): Breath rate, breath waveform; 6) invasive blood pressure (IBP): 2 channels IBP waves and Systolic pressure, diastolic pressure, mean pressure;
- carbon dioxide (CO₂): 1 channel CO₂ wave and Breath rate, Et CO₂ and FiCO₂.

Features

- 13 Type arrhythmic analysis, multi lead ECG waveforms display in phase, real time ST segment analysis, pacemaker detection drug calculation and titration table
- Efficient resistance to interference of defibrillator and electrosurgical cautery;
- SPO₂ can testing for 0.1% weak;
- RA-LL impedance respiration;
- Trend coexist display;
- Options of printing, IBP, VGA and EtCO₂;
- UP to 4 hours working capacity of built_in rechargeable battery;
- 15" high resolution color TFT LCD display;
- Large volume of tabular and graphic trends information storage and easy to recall ;
- Anti-ESU, anti- defibrillator;
- Capture dynamic waveforms

- OxyCRG DyNamic view display