


Clipper



Clipper is an endoscopic instrument for imposing titanium clips with a length of 3 mm to 11 mm on anatomical tubular elements: vessels, ducts, etc.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Clipper is an endoscopic instrument for imposing titanium clips with a length of 3 mm to 11 mm on anatomical tubular elements: vessels, ducts, etc.

NO.	Name	Characteristics
CA0001L	Clip applicator, medium with one and two movable jaws	Diameter 10 mm, working length 330 mm
CA0002L	Clip applicator, small	Diameter 5 mm, working length 330 mm
CA0003L	Clip applicator, medium	Diameter 10 mm, working length 330 mm
CA0004L	Clip applicator for deep layers	Diameter 10 mm, working length 330 mm
CA0005L	Clip applicator for plastic clips, medium	Diameter 5 mm, working length 330 mm
CA0006L	Clip applicator for plastic clips, large	Diameter 10 mm, working length 330 mm
CA0007L	Clip applicator for plastic clips, small	Diameter 5 mm, working length 330 mm
CA0008L	Clip Applicator	
CA0009L	Clip applicator for biological clips	
TC0001L	Titanium clips for laparoscopy	
PC0001L	Plastic clips for laparoscopy	