

## Diagnostic set LTI-100



Diagnostic set large

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

Diagnostic set large (Bayonet connection system)

Lightweight, robust zipper case containing one 2.5 v LED Parker otoscope, 3 ear speculum (size 2.5 3.5 4.5mm), 1 Parker ophthalmoscope, 1 bend arm throat lamp, 2 laryngeals mirrors sizes 3 & 4, 1 tongue spatula, support for single use tongue depressor, 1 adjustable nasal speculum fitting, 2 bottles spare mini specula, black cap bottle have 10 specula size 2.5mm and white cap bottle have 10 specula size 4mm and 1 spare bulb 2.5v mm for otoscope.