

Diagnostic set LTI-105



Diagnostic set (Bayonet connection system)

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Diagnostic set (Bayonet connection system).

- Otoscope head: pin contacting fitting.
- Chrome plated fixed lens 4X in plastic mount may be removed for introductions for instruments.
- Reusable specula: (2.5 3.5 4.5mm) included.
- Ophthalmoscope head: pin contact fitting.
- Lens selection wheel permits choice of corrective lenses from 0 to -25 and 0 to +40 diopters single large circle aperture.
- Bent-arm illumination
- Tongue blade holder
- Plastic Nasal speculum
- Chrome plated
- Two laryngeal mirrors: #3 and #4
- Handle Chrome plated
- Rheostat for adjustment of illumination