

## DSS Xenon Light Source



The Brightfield Model DSS01 Xenon Light Source is designed to work with endoscopes, imaging systems, and other ancillary equipment for imaging-guided diagnostic and therapeutic endoscopic procedures.

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

The Brightfield DSS High Intensity Xenon Light Source is designed for beam shaping and illumination of the surgical site for endoscopic procedures in a variety of clinical areas.

The light source in the illuminator is a xenon lamp, which has a radiation spectrum as close as possible to natural. The lamp resource is more than 500 hours. The DSS model illuminator allows you to obtain exceptional illumination of the affected area with a greater efficiency.

### Features of DSS Xenon Light Source

- A modern xenon lamp with stable temperature characteristics and high brightness.
- The pulse width rectifier provides power at the specified output data.
- Good accessibility for changing the light source is provided.
- Illumination field of view 1200000 lx
- Color temperature 2 300 K
- Weight 4.5 kg