

Endolaryngeal surgery instruments



Standard Instruments for Transoral Endolaryngeal Surgery

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

Standard Instruments for Transoral Endolaryngeal Surgery.

NO.	Product name	Characteristics
BF0001L	Biopsy forceps with spoon-shaped jaws	Diameter 2 mm, long 230 mm
BF0002L	Biopsy forceps with spoon-shaped jaws curved upwards at 45 degrees	Diameter 2 mm, long 230 mm
BF0003L	Biopsy forceps with spoon-shaped jaws curved to the left	Diameter 2 mm, long 230 mm
BF0004L	Biopsy forceps with spoon-shaped jaws curved to the right	Diameter 2 mm, long 230 mm
BF0005L	Biopsy forceps with spoon-shaped elongated jaws	Diameter 2 mm, long 230 mm
BF0006L	Biopsy forceps with spoon-shaped jaws, bent upwards by 40 degrees	Diameter 2 mm, long 230 mm
BF0007L	Biopsy forceps with spoon-shaped jaws bent to the left by 45 degrees	Diameter 2 mm, long 230 mm
BF0008L	Biopsy forceps with flat spoon-shaped jaws bent to the right at 45 degrees	Diameter 2 mm, long 230 mm
BF0009L	Biopsy forceps with flat spoon-shaped jaws, bent upwards by 40 degrees, large	Diameter 3mm, long 230mm
BF0010L	Biopsy forceps with flat spoon-shaped jaws, curved upwards by 40 degrees, medium	Diameter 3 mm, long 230 mm
BF0011L	Biopsy forceps with flat spoon-shaped jaws, curved upwards by 40 degrees, small	Diameter 2 mm, long 230 mm
FF0001L	Foreign body pliers, denticulate	Diameter 2 mm, long 230 mm
SM0001L	Miniature straight scissors.	Diameter 2 mm, long 230 mm