

Endoscopic laryngoscope



An endoscopic laryngoscope is a rigid endoscopic instrument designed to examine the pharyngeal cavity and vocal cords and conduct minimally invasive interventions in this anatomical area.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

An endoscopic laryngoscope is a rigid endoscopic instrument designed to examine the pharyngeal cavity and vocal cords and conduct minimally invasive interventions in this anatomical area.

Features of Brightfield Rigid Laryngoscopes

- The body of the tube and its constituent elements are made of high grade stainless steel
- The endoscope uses elements of optics from a German manufacturer.
- The optical axes are created using the patented lens shaft technology, which ensures the visibility of the entire field of view.
- The cap and sapphire lens are durable and never wear off.
- The valve is made of stainless steel, easy to maintain, protected from damage.
- The connectors are compatible with endoscopic equipment from leading manufacturers, including Storz products.
- Sterilization can take place at low temperatures.

Brightfield laryngoscope optical systems

NO.	Field of view	Diameter	Working length	Note
LA_01	0 degrees	10 mm	310 mm	
LA_02	0 degrees	10 mm	325 mm	
LA_03	0 degrees	10 mm	450 mm	
LA_04	0 degrees	10 mm	540 mm	
LA_05	0 degrees	11 mm	274 mm	angled eyepiece
LA_06	0 degrees	5.5 mm	310 mm	
LA_07	0 degrees	5 mm	290 mm	angled eyepiece
LA_08	0 degrees	5 mm	175 mm	
LA_09	0 degrees	5 mm	300 mm	angled eyepiece
LA_10	0 degrees	5 mm	300 mm	
LA_11	0 degrees	5 mm	360 mm	
LA_12	0 degrees	6.5 mm	300 mm	
LA_13	0 degrees	4 mm	185 mm	angled eyepiece
LA_14	30 degrees	10 mm	310 mm	
LA_15	30 degrees	10 mm	325 mm	

Endoscopes: Endoscopic laryngoscope

LA_16	30 degrees	10 mm	450 mm	
LA_17	30 degrees	10 mm	540 mm	
LA_18	30 degrees	5.5 mm	310 mm	
LA_19	30 degrees	5 mm	240 mm	
LA_20	30 degrees	5 mm	290 mm	
LA_21	30 degrees	5 mm	300 mm	
LA_22	30 degrees	5 mm	360 mm	
LA_23	30 degrees	6.5 mm	300 mm	
LA_24	45 degrees,	10 mm	310 mm	
LA_25	45 degrees	10 mm	325 mm	
LA_26	45 degrees	10 mm	450 mm	
LA_27	45 degrees	5 mm	300 mm	
LA_28	12 degrees	4.8 mm	183 mm	angled eyepiece
LA_29	12 degrees	4 mm	185 mm	
LA_30	70 degrees	10 mm	186 mm	
LA_31	70 degrees	11 mm	200 mm	
LA_32	70 degrees	4 mm	175 mm	
LA_33	70 degrees	5 mm	175 mm	
LA_34	70 degrees	5 mm	240 mm	
LA_35	70 degrees	6 mm	186 mm	with a handle
LA_36	70 degrees	6 mm	186 mm	
LA_37	70 degrees	6 mm	190 mm	
LA_38	70 degrees	8 mm	170 mm	
LA_39	70 degrees	8 mm	186 mm	with a handle
LA_40	70 degrees	8 mm	186 mm	
LA_41	70 degrees	8 mm	190 mm	
LA_42	90 degrees	10 mm	190 mm	
LA_43	90 degrees	11 mm	200 mm	
LA_44	90 degrees	6 mm	180 mm	
LA_45	90 degrees	6 mm	186 mm	with a handle
LA_46	90 degrees	6 mm	186 mm	
LA_47	90 degrees	8 mm	174 mm	
LA_48	90 degrees	8 mm	178 mm	
LA_49	90 degrees	8 mm	180 mm	
LA_50	90 degrees	8 mm	186 mm	with a handle
LA_51	90 degrees	8 mm	186 mm	
LA_52	90 degrees	8 mm	190 mm	
LA_53	12 degrees	4.8 mm	183 mm	laryngoscope with working channel and angled eyepiece