

Forceps and scissors for proctology.



For the collection of material for histological studies and minimally invasive interventions in the rectal area, most of the use is made of flexible or rigid grasping serrated forceps with two movable jaws, as well as dissection (excisional) forceps of various configurations.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

For the collection of material for histological studies and minimally invasive interventions in the rectal area, most of the use is made of flexible or rigid grasping serrated forceps with two movable jaws, as well as dissection (excisional) forceps of various configurations.

NO.	Name	Characteristics
BC0001FL	Flexible biopsy forceps with two movable jaws	diameter 7 Fr, working length 330 mm
GC0001FL	Flexible serrated pliers with two movable jaws	diameter 7 Fr, working length 330 mm
EC0001FL	Flexible toothed foreign body forceps with two movable jaws	diameter 7 Fr, working length 330 mm
SC0001FL	Flexible scissors with one movable blade	diameter 7 Fr, working length 330 mm
CP0001P	Flexible pliers with two movable jaws	