

Gastroscope



The gastroscope is a flexible endoscopic instrument for diagnostic and therapeutic manipulations in the esophagus, stomach and duodenum 12.

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

The gastroscope is a flexible endoscopic instrument for diagnostic and therapeutic manipulations in the esophagus, stomach and duodenum 12. Today on the market of endoscopic equipment there are gastroscopes of several designs - a video gastroscope and a fibrogastroscope. Visualization on a video gastroscope is due to a miniature digital matrix. The image during fibrogastroscopy is transferred through the fiberglass directly to the monitor screen.

Features of Brightfield gastroscopes

- Superior image quality and fidelity
- Visualization is carried out using a modern digital matrix
- Increased field of view
- Optimal range of motion for the distal end
- Compatible with endoscopic instruments of all manufacturers
- Connectors are compatible with endoscopic equipment from leading manufacturers

NO.	Field of view	Depth of field	Inserted tube diameter	Working channel diameter	total length
GAST_FL04	120 degrees	3 - 100 mm	9.2 mm	2.8 mm	1300 mm