

Handles



Handles for laparoscopic instruments with a diameter of 5 and 10 mm. Tool assembly is fast and reliable. The use of collapsible tools greatly facilitates the cleaning process and extends the service life of the tool.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Handles for laparoscopic instruments with a diameter of 5 and 10 mm. Tool assembly is fast and reliable. The use of collapsible tools greatly facilitates the cleaning process and extends the service life of the tool.

NO.	Name	Characteristics
HI0005L	Handle for laparoscopic instruments	Diameter 5 mm
HI0010L	Handle for laparoscopic instruments	Diameter 10 mm
HI0S10L	Handle for laparoscopic instruments with lock	Diameter 10 mm
NH0001L	Ring handle for needle holders	Diameter 5 mm
NH0002L	Straight, V-shaped handle for needle holders	Diameter 5 mm