

## ICU Ventilator S1400



The ICU Ventilator S1400 with an 18.5-inch screen and innovative ventilation modes PC-APRV and PC-ACV. The device is designed to work in all types of medical facilities where there is a need for forced ventilation.

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

The ICU Ventilator S1400 with an 18.5-inch screen and innovative ventilation modes PC-APRV and PC-ACV. The device is designed to work in all types of medical facilities where there is a need for forced ventilation.

Ventilation modes: MANU, SPONT, VC-ACV, VC-SIMV, PSV, SIMV PCV, Spn-CPAP, PCV, PC-APRV, PC-ACV, PRVC, S-PRVC

### Features

- 18.5-inch touch screen large screen with a clear display of all functions;
- Display of 5 curves and 3 loops;
- 12 types of ventilation modes;
- Three-level alarm system, visual and audible alarm;
- Self-test before starting, eliminating a system error;
- Improved power system control technology;
- Availability of built-in backup power supply; when the external power source is discharged, the backup power source turns on automatically.
- Separate placement of electronic module and gas supply system.