

## Instruments for proctoscopy



A set of reusable instruments for examining the lower rectum and anus using a fiber-optic light source. The working element of the proctoscope has an instrumental channel for minimally invasive visual control interventions

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

A set of reusable instruments for examining the lower rectum and anus using a fiber-optic light source. The working element of the proctoscope has an instrumental channel for minimally invasive visual control interventions.

NO.	Name	Characteristics
OB0001P	Obturator with tube	diameter 20 mm, working length 123 mm
RP0001P	Work item	
CG0001P	Glass cap	
ST0001P	Curved suction tube	
ST0001P	Cleaning rod	
IN0001P	Insufflation bulb with connector	