


Knot pusher



Endoscopic instrument for the stitching and tying of various sutures in laporoscopic operations

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Endoscopic instrument for the stitching and tying of various sutures in laporoscopic operations.

NO.	Name	Characteristics
KF0001L	Pulling and knotting pliers, rotary	Diameter 5 mm, working length 330 mm
KF0002L	Pulling and knotting pliers, type B Diameter 5 mm, working length 100 mm	Diameter 5 mm, working length 330 mm
KP0001L	Knot pusher	Diameter 5 mm, working length 330 mm