

Laryngoscope LTI-7000



The Flex tip is a Macintosh laryngoscope blade

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

The Flex tip is a Macintosh laryngoscope blade

The Flex tip is a Macintosh laryngoscope blade with a tip which is adjustable through 70%. A lever controls the tip angle during incubation to lift the epiglottis giving a clear view of the vocal cords. Contact with upper row of teeth and unnecessary pressure on the tongue and soft tissue and avoided.

- Swivelling blade tip facilitates difficult incubation.
- Reduced likelihood of damage to teeth and soft tissue.
- Integrated Fiber Optics for the best possible illumination.
- High quality stainless steel construction.
- Fits all standard EO. handles (green standard)