


Laser surgery



A nephroscope is an endoscopic instrument that is designed for examination and manipulation of the renal pelvis

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Optical systems for the treatment of urolithiasis and soft tissue therapy using thulium and holmium lasers.

Laser therapy techniques are relatively simple, do not require expensive equipment, and are effectively combined with almost all other methods of treatment, both therapeutic and surgical. This explains the widespread use of the technique in modern clinical practice.

Currently, more than 200 European centers have laser equipment for one prostatectomy and other soft tissue interventions. New methods and technologies continue to appear in urology, which take advantage of the advantages of the laser over traditional treatment.

NO.	Field of view	Diameter	Length
ULAS01	12 degrees	4 mm	302 mm
ULAS02	30 degrees	4 mm	302 mm
REZ01	12 degrees	4 mm	302 mm