

Morcellator



A set of instruments for dissection and removal of large organs from the connective muscle tissue during gynecological laparoscopic interventions, such as excision of fibroids, supracervical hysterectomy, etc.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

A set of instruments for dissection and removal of large organs from the connective muscle tissue during gynecological laparoscopic interventions, such as excision of fibroids, supracervical hysterectomy, etc.

NO.	Name	Characteristics
ME1	The electronic unit	
ME2	Electro-mechanical unit	
KN0001M	Morcellator knife	Diameter 10 mm, working length 320 mm
KN0002M	Morceller knife	Diameter 15 mm, working length 320 mm
KN0003M	Morceller knife	Diameter 18 mm, working length 320 mm
KN0004M	Morceller knife	Diameter 15 mm, working length 320 mm
KN0005M	Morceller knife	Diameter 18 mm, working length 320 mm
ZN0001M	Measuring probe, graduated	Diameter 5 mm, working length 500 mm
ZN0002M	Probe	Diameter 5 mm, working length 420 mm.
ZN0003M	Probe conductor	Diameter 10 mm, working length 330 mm.
KF0001M	Knot pusher	Diameter 5 mm, working length 330 mm.
HB001M	Myoma fixator	Diameter 5 mm, working length 420 mm.
HB002M	Myoma fixator	Diameter 10 mm, working length 420 mm.
HB003M	Myoma fixator	Diameter 5 mm, working length 330 mm.
DL002M	Dilator	Diameter 15 mm
DL001M	Dilator	Diameter 18 mm
CO001M	Converter	Diameter 15 - 10 mm
CO002M	Converter	Diameter 18 - 10 mm