

Operating lights LED160+80

Double operating lights 160 LED and 80 LED

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Double operating lights 160 LED and 80 LED

- Lamp diameter, mm: 830 + 640
- Lamp type: convex lens
- Maximum illumination, Lux: 180000 + 120 000
- Colour temperature, K: 4300
- CRI: more 95
- Light field diameter, mm: 100 - 300
- Light field depth, mm: 850
- Number of LEDs: 160 EA + 80 EA
- LED life time: 50 000 hour
- Touch on/off: touch panel & remote control (option)
- Power consumption, W: 240 + 120
- Dimming range: 20-100% (5 steps)
- Opto semiconductors: Germany medical Osram
- Spring arm: Germany Ondal

