

Operating lights LED56+40

Double operating lights 56 LED and 40 LED

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Double operating lights 56 LED and 40 LED

- Lamp diameter, mm: 685 + 550
- Lamp type: convex lens
- Maximum illumination, Lux: 160000 + 120 000
- Colour temperature, K: 4300
- CRI: more 95
- Light field diameter, mm: 100 - 300
- Light field depth, mm: 850
- Number of LEDs: 96 (56 + 40) EA
- LED life time: 50 000 hour
- Touch on/off: touch panel & remote control (option)
- Power consumption, W: 84 + 60
- Dimming range: 20-100% (5 steps)
- Opto semiconductors: Germany medical Osram
- Spring arm: Germany Ondal

