

Otoscope LTI-2009



High quality professional otoscope

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

High quality professional otoscope

High quality professional plastic head, chrome plated battery handle Otoscope with extra bright light 2.5v LED bulb 2.5v illumination provided through fiber optic system. High quality big swiveling viewing window with 3x magnification lens out put for pneumatic membrane test.

Supply: 10 single use mini speculum (0 2.5 and 4mm) with storage bottle packing in black pouch.

Require: 2x C batteries (not included with set)