

## Endoscopy Pump - PERISTALTIC LAP



Endoscopy Pump PERISTALTIC - for Laparoscopy

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

Endoscopy Pump PERISTALTIC - for Laparoscopy

- Separate channel for suction and for irrigation
- Irrigation through the peristaltic pump module
- Suction through the vacuum pump module
- Predefined suction and irrigation parameters

**Pressure range:** 350-450 mmHg

**Flow range:** 1,8 l/min

**Suction pressure:** -50 kPa to -60 kPa

**Dimensions (DxWxH):** 310 x 330 x 120 mm

**Weight:** 8 kg

**Supply voltage:** 100-240 V AC / 50/60 Hz

**Class of protection:** Applied part - BF type

**Casing safety standart:** IP X1

**CE Standart:** N60601-1;EN60601-1-2