

## GL350-5A dual lamp light source



The Brightfield GL350-5A high intensity light source is equipped with a xenon and halogen lamp.

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

The Brightfield GL350-5A high intensity light source is equipped with a xenon and halogen lamp. Such a design solution allows the device to be used both in laparoscopy, where a powerful light source is required, and in thoracic endoscopic operations, where it is possible to limit ourselves to a halogen light source.

### Features light source model GL350-5A

- Xenon lamp power 350 W
- Halogen lamp power 250
- Color temperature 5500 K
- Power supply: 220V
- Frequency 50-60 MHz