

Laryngoscope



Laryngoscope is a rigid endoscopic instrument for diagnostic studies and minimally invasive interventions in the upper respiratory tract (pharynx, larynx, vocal cords, etc.)

Rating: Not Rated Yet
[Ask a question about this product](#)

Description

Laryngoscope is a rigid endoscopic instrument for diagnostic studies and minimally invasive interventions in the upper respiratory tract (pharynx, larynx, vocal cords, etc.). Brightfield offers a wide range of endoscopic laryngoscopes with various optical systems and instruments for surgical intervention.

Features of Brightfield endoscopic laryngoscopes

- The body of the tube and its constituent elements are made of high grade stainless steel
- The endoscope uses elements of optics from a German manufacturer.
- The optical axes are created using the patented lens shaft technology, which ensures the visibility of the entire field of view.
- The cap and sapphire lens are durable and never wear off.
- The valve is made of stainless steel, easy to maintain, protected from damage.
- The connectors are compatible with endoscopic equipment from leading manufacturers, including Storz products.
- Sterilization can take place at low temperatures.

NO.	Field of view	Diameter	Working length	Note
LATCH_01	0 degrees	4 mm	185 mm	angled eyepiece
LATCH_02	12 degrees	4.8 mm	183 mm	angled eyepiece
LATCH_03	12 degrees	4 mm	185 mm	angled eyepiece
LATCH_04	30 degrees	5 mm	240 mm	
LATCH_05	70 degrees	10 mm	186 mm	
LATCH_06	70 degrees	11 mm	200 mm	
LATCH_07	70 degrees	4 mm	175 mm	
LATCH_08	70 degrees	5.8 mm	175 mm	
LATCH_09	70 degrees	5 mm	175 mm	
LATCH_10	70 degrees	5 mm	240 mm	
LATCH_11	70 degrees	6 mm	186 mm	with a handle
LATCH_12	70 degrees	6 mm	186 mm	
LATCH_13	70 degrees	6 mm	190 mm	
LATCH_14	70 degrees	8 mm	170 mm	
LATCH_15	70 degrees	8 mm	186 mm	with a handle
LATCH_16	70 degrees	8 mm	186 mm	

LATCH_17	70 degrees	8 mm	190 mm	
LATCH_18	90 degrees	10 mm	190 mm	
LATCH_19	90 degrees	11 mm	200 mm	
LATCH_20	90 degrees	6 mm	180 mm	
LATCH_21	90 degrees	6 mm	186 mm	with a handle
LATCH_22	90 degrees	6 mm	186 mm	
LATCH_23	90 degrees	8 mm	174 mm	
LATCH_24	90 degrees	8 mm	178 mm	
LATCH_25	90 degrees	8 mm	180 mm	with a handle
LATCH_26	90 degrees	8 mm	186 mm	with a handle
LATCH_27	90 degrees	8 mm	186 mm	