

Vein Illuminator BR-1.01



Medical vein illuminator for hospital

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Medical vein illuminator for hospital

We offer hospitals the most innovative peripheral vein illumination products on the market. With HD imaging and exclusive DLP technology, our vein finder assists your hospital staff in seeing deeper vasculature with more accuracy than other vein light devices on the market.

Basing on the principle that hemoglobin in the blood vessel has stronger absorption of near-infrared light than other tissues and through a series of digital image processing, Projection Vein Finder projects the outline of the vein to the skin surface. It provides a real time accurate image of the patient's blood pattern so that medical personnel can clearly observe the tiny blood vessels under the skin. It greatly reduces the difficulty of venous puncture and patients' fear of injections.