

Trocars



Endoscopic instrument for penetration into the abdominal cavity through the integument with the preservation of the instrumental canal and its sealing.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Endoscopic instrument for penetration into the abdominal cavity through the integument with the preservation of the instrumental canal and its sealing.

NO. Name	Characteristics
TR Magnetic Valve 000 Trocar 1L	Diameter 5.5 mm, working length 100 mm
TR Magnetic Valve 000 Trocar 2L	Diameter 10.5 mm, working length 100 mm
TR Protective solenoid 000 trocar 3L	Diameter 5.5 mm, working length 100 mm
TR Protective solenoid 000 trocar 4L	Diameter 10.5 mm, working length 100 mm
TR Screw thread 000 solenoid trocar 5L	Diameter 5.5 mm, working length 100 mm
TR Screw thread 000 solenoid trocar 6L	Diameter 10.5 mm, working length 100 mm
TR Trocar, cruciform 000 7L	Diameter 5.5 mm, working length 100 mm
TR Trocar, cruciform 000 8L	Diameter 10.5 mm, working length 100 mm
TR Trocar, cruciform 000 with protection 9L mechanism	Diameter 5.5 mm, working length 100 mm

Tool Part: Trocars

TR Trocar, cruciform 001 with protection 0L mechanism	Diameter 5.5 mm, working length 100 mm
TR Trocar, cruciform 001 with screw thread 1L	Diameter 5.5 mm, working length 100 mm
TR Trocar, cruciform 001 with screw thread 2L	Diameter 10.5 mm, working length 100 mm
TR Trocar, cruciform, 001 variable 3L	Diameter 10.5 mm, working length 100 mm
TR Magnetic Valve 001 Trocar with Vent 4L Valve	Diameter 5.5 mm, working length 100 mm
TR Magnetic Valve 001 Trocar with Vent 5L Valve	Diameter 10.5 mm, working length 100 mm
TR Magnetic Valve 001 Trocar with Vent 6L Valve with Protection Mechanism	Diameter 5.5 mm, working length 100 mm
TR Magnetic Valve 001 Trocar with Vent 7L Valve with Protection Mechanism	Diameter 10.5 mm, working length 100 mm
TR Solenoid trocar with 001 a vent valve with 8L screw thread	Diameter 10.5 mm, working length 100 mm
TR Variable solenoid 001 trocar with vent 9L valve	Diameter 10.5 mm, working length 100 mm
TR Vent Valve Trocar 002 0L	Diameter 5.5 mm, working length 100 mm
TR Vent Valve Trocar 002 1L	Diameter 10.5 mm, working length 100 mm
TR Trocar with 002 ventilation valve with 2L protection mechanism	Diameter 5.5 mm, working length 100 mm
TR Trocar with 002 ventilation valve with 3L protection mechanism	Diameter 10.5 mm, working length 100 mm
TR Threaded trocar 002 4L	Diameter 12 mm
TR Threaded trocar 002 5L	Diameter 10 mm
TR Threaded trocar 002 6L	Diameter 5 mm
RA Transition cannula 000 1L	Diameter 10.5 Fr to 5.5 Fr
RA Syringe tip 000 2L	Diameter 10.5 / 5.5 mm
NV Veress needle 000 1L	Diameter 2.2 mm, working length 120 mm