


Tubes and seals



External and internal tubes, obturators and other elements of the endoscopic system for examining the vagina, cervix and uterine cavity

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

External and internal tubes, obturators and other elements of the endoscopic system for examining the vagina, cervix and uterine cavity

NO.	Product name	Specifications
TI0001G	Inner tube & Obturator	12 Fr. x 232 mm
TI0002G	Inner tube & Obturator	16 Fr. x 232 mm
TO0003G	Outer Sheath & Obturator	25 Fr x 195 mm
TO0004G	Outer Sheath & Obturator	16 Fr x 209 mm
OB0005G	Obturator	12 Fr
TG0001G	Tube	19 Fr. x 208 mm
TG0002G	Tube	15 Fr. x 209 mm
TI0003G	Tube	15 Fr. x 239 mm
TG0004G	Tube	27 Fr x 180 mm
SW0001G	Steel conductor	0.8 Fr. x 300 mm
RC0001G	Cleaning rod	6.0 Fr.
WG0001G	Work item	19.5 Fr x 190 mm
WG0002G	Work item	16 Fr x 198 mm
WG0003G	Work item	14 Fr x 191.5 mm
WD0004G	Working element with external protection	13 Fr x 220 mm
WP0005G	Work item (passive)	
WA0006G	Work item (active)	
VI0001	Inlet flow valve	
VO0001	Outlet flow valve	
VH0001	Hermetic cap	