

UProbe-2L



The new generation instruments for ultrasonography, redefining access in the interventional suite and at the point of care

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

The new generation instruments for ultrasonography, redefining access in the interventional suite and at the point of care

UProbe-2L - the new generation instruments for ultrasonography, redefining access in the interventional suite and at the point of care. The probe acted as a Wi-Fi Access Point and can be connected by iPad. With the probe be connected through WiFi and the App is running, enjoy your days of working without the trouble making cables.

Application: thyroid, breast, vascular, nerve, MSK

Specification:

- Scan mode: electronic array
- Length: 40 mm
- 64 elements, 8 channel
- Display mode: B, B/M
- Frequency?7.5MHz
- Depth: 100 – 200 mm
- Measure: distance, area, obstetrics
- Battery last: 4 hours
- Dimension: 180mmx60mmx60mm
- Weight: 300 gram
- Wi-Fi type: 802.11g/20MHz/2.4G/54Mbps
- Work OS: Apple iOS and Android
- Image frame rate: 12f/s