

UProbe-3C



Wireless convex ultrasound transducer - handheld wireless ultrasound machine.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Wireless convex ultrasound transducer - handheld wireless ultrasound machine.

Wireless convex ultrasound transducer UProbe-3C - handheld wireless ultrasound machine. It offers high resolution ultrasound images that are as good as the best traditional point-of-care systems.

- Workable with Tablet or Smart Phone
- Built-in and replaceable battery
- Advanced digital imaging technology, clear image
- Wireless connectivity, easy to operate
- Small and light, easy to carry
- Intelligent terminal platform, powerful expansion functions on application, storage, communication, printing
- Application: abdomen, obstetric

Specification:

- Scan mode: electronic array
- Radius: 60 mm
- 128 elements, 16 channel
- Display mode: B, B/M
- Frequency?3.5 MHz
- Depth: 100 – 200 mm
- Measure: distance, area, obstetrics
- Battery last: 3 hours
- Dimension: 156mm×60mm×20mm
- Weight: 308 gram
- Wi-Fi type: 802.11g/20MHz/2.4G
- Work OS: Apple iOS and Android
- Image frame rate: 12f/s