

UProbe-3N



This wireless probe scanner with integrated ultrasound circuit boards inside

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

This wireless probe scanner with integrated ultrasound circuit boards inside

This wireless probe scanner UProbe-3N with integrated ultrasound circuit boards inside, can wireless connect to the PC Tablet, Smartphone that after installed ultrasound software, then realize the functions of a ultrasound scanner. It is so small and smart, convenient for carry and operation.

Application: in emergency clinical, hospital ward inspection, community clinical and outdoor inspect

Specification:

- Scan mode: electronic array
- Length: 25 mm
- 128 elements, 16 channel
- Display mode: B, B/M
- Frequency?10 MHz
- Depth: 100 – 200 mm
- Measure: distance, area, obstetrics
- Battery last: 4 hours
- Dimension: 180mmx60mmx60mm
- Weight: 300 gram
- Wi-Fi type: 802.11g/20MHz/2.4G/54Mbps
- Work OS: Apple iOS and Android
- Image frame rate: 12f/s