

## UProbe-4L



The new generation instruments for ultrasonography with the outstanding feature of wireless

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

The new generation instruments for ultrasonography with the outstanding feature of wireless

UProbe-4L - the new generation instruments for ultrasonography with the outstanding feature of wireless. The probe is highly integrated with ultrasound image processing, power management and a wireless signal provider to be connected by the main unit.

The main units different with traditional devices are now changed to be any iPad from Apple Inc. The probe acted as a Wi-Fi Access Point and can be connected by iPad. With the probe be connected through WiFi and the App is running, enjoy your days of working without the trouble making cables.

Application: thyroid, breast, vascular, nerve, MSK

### Specification:

- Scan mode: electronic array
- Length: 40 mm
- 128 elements, 32 channel
- Display mode: B, B/M
- Frequency?7.5MHz
- Depth: 100 – 200 mm
- Measure: distance, area, obstetrics
- Battery last: 4 hours
- Dimension: 180mmx60mmx60mm
- Weight: 300 gram
- Wi-Fi type: 802.11g/20MHz/2.4G/54Mbps
- Work OS: Apple iOS and Android
- Image frame rate: 12f/s