

Vein Illuminator BR-1.02



Medical vein illuminator for hospital

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

New integrated exam light

Basing on the principle that hemoglobin in the blood vessel has stronger absorption of near-infrared light than other tissues and through a series of digital image processing, Projection Vein Finder projects the outline of the vein to the skin surface, so that medical personnel can clearly observe the tiny blood vessels 8mm to 10mm under the skin.

- Adjustable Brightness: Through the operation panel, the brightness of projection images can be adjusted for different light sensitivities of personnel.
- Adult& Child Mode: Projection image area is different according to different groups, It is user-friendly.
- Portable and Easy to Adjust: B500 has a size of 227mm*65mm*63mm (length * width * height) and weight about 500g, It is portable in the pocket.
- Clamp Support turns around in 360°, and adjusts up and down.
- It has function of vein depth recognition.