

Vein Illuminator (pediatric)



The battery-powered pediatric medical vein illuminator

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

The battery-powered pediatric medical vein illuminator

The battery-powered pediatric medical vein illuminator emits a special wavelength of red light, which can penetrate the patient's extremities, and after absorbed by blood, the veins become darker than the surrounding tissues. There for the contrast ratio between the veins and the surrounding tissues is increased. It assists the paramedics to easily locate the vein and observe extravasated blood for patients' extremities. To prevent cross-infection in clinic between patients or doctors and patients, this product provides a special isolation cover.

- Lighting Power < 20 VA
- Battery Duration > 2h
- Power Supply 4V

