

Veinlite



Medical vein illuminator 28 bright light emitting diodes

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Medical vein illuminator 28 bright light emitting diodes

The adult model with optimum viewing area and highest contrast for vein imaging. The rechargeable Veinlite features 28 bright light emitting diodes (LEDs). The 22 bright orange and 6 deep red LEDs can be operated together or separately for maximum versatility and use on patients of all ages and skin tones. With this model there is no need to keep spare batteries in hand.

- Fewer IV access failures
- Fewer needle sticks
- Less pain in for the patient
- Fewer wasted supplies